

PRECURSOR ION INTENSITY OR SPECTRAL COUNTING?

LC-MS/MS data processing and statistical assessment for relative protein-level quantification in label-free proteomics

ANANDASHANKAR ANIL, ANGELIKI MARAKI, JAMES ERICSON

aanil@kth.se, aggmara@gmail.com, jameser@kth.se

Precursor Peak Intensity Method

- For each sample:
 - ❖ Normalization of the Intensity Area of each peptide against the total Intensity of the sample.
 - ❖ Sum normalized intensity for all identical peptides identifying a protein.
 - Unique peptides (summed intensity area) identifying a single protein
 - Removal of peptides that do not occur in all three replicates
- Comparison between two samples:
 - ❖ Create a single list containing all peptides that occur in both samples
 - ❖ Perform **t.test** on $\log(\text{sum}(\text{intensities}))$
 - ❖ Sort ascending by p-value from t.test
 - ❖ Obtain **q-value**

Significantly Differentially Abundant Proteins

1v2	
Protein Name	
sp P55249 ZRT4_YEAST	
sp P38788 SSZ1_YEAST	
sp P44983 UTR6_YEAST	

2v4	
Protein Name	
sp P44983 UTR6_YEAST	
sp P00927 THDH_YEAST	
sp P44374 SFG2_YEAST	
sp P40825 SYA_YEAST	

3v4	
Protein Name	
sp P44683 PGA4_YEAST	
sp P32473 ODPB_YEAST	
sp P44983 UTR6_YEAST	
sp P40825 SYA_YEAST	
sp Q05016 YM71_YEAST	
sp P19097 FAS2_YEAST	

1v4	
Protein Name	
sp P44683 PGA4_YEAST	
sp P44374 SFG2_YEAST	
sp P40825 SYA_YEAST	
sp P44983 UTR6_YEAST	
sp Q05016 YM71_YEAST	
sp P11075 SEC7_YEAST	

Significantly Differentially Abundant Proteins

1v3	
Protein Name	
sp P44374 SFG2_YEAST	

2v3	
Protein Name	
sp P55249 ZRT4_YEAST	

Spectral Counting Method

- NSAF(Normalized Spectral Abundance Factor)
 - ❖ Filter: get only the peptide measurements for which the q value was less than 0.01.
 - ❖ NSAF calculated as the number of spectral counts(SpC) identifying a protein, divided by the protein's length(L), divided by the sum of SpC/L for all proteins in the replicate
 - ❖ Merged the NSAF scores for all replicates in two specific runs and performed a t.test.
 - ❖ Obtained q-value
- emPAI(Exponentially modified protein abundance index)
 - ❖ Filter: get only the peptide measurements for which the q value was less than 0.01.
 - ❖ Removed [0-9] modification sites in peptides for comparison against digested peptides
 - ❖ Removed peptide duplicates and count Observed peptides
 - ❖ Calculated normalized emPAI by dividing individual emPAI score with sum of all emPAI scores from a replicate in a sample
 - ❖ Merged the emPAI scores for all replicates in two specific runs and performed a t.test.
 - ❖ Obtained q-value

No Significantly Differentially Abundant Proteins Found

Conclusion

- ❖ The Precursor ion peak intensity method identified more differentially abundant proteins than the Spectral Counting methods.
- ❖ Out of the two Spectral Counting approaches, the NSAF method identified only two significantly differentially abundant proteins, while emPAI none.